

## Thermal expansion and magnetostriction measurements

An application using Oxford Instruments Superconductivity's HelioxVL, advanced sorption single shot  $^3\text{He}$  refrigerator

The samples available for such measurements are frequently of the order of  $1\text{ mm}^3$  only. By using tilted capacitor plates the new sensor can measure samples with an active length of less than  $1\text{ mm}$ , with an absolute resolution of about  $1\text{ \AA}$ . The dilatometer (Figure 1) is designed to fit into a variable temperature insert with a  $50\text{ mm}$  bore or less. For low temperature measurements between  $300\text{ mK}$  and  $4\text{ K}$  an Oxford Instruments Superconductivity HelioxVL  $^3\text{He}$  insert is used.

Dr Rotter and his colleagues (Figure 2) used the tilted plate method to overcome some of the problems associated with parallel plates. The capacitance method is one of the most sensitive methods for measuring small length changes of solids and has the advantage, compared with specific heat measurements, of covering a large temperature range with one instrument. Capacitive cells with parallel plates are easy to calibrate, but they are either large, causing problems with thermal stability, or present difficulties with sample handling. Even if the problem of thermal stability can be overcome, the limited space available in the magnetic coil causes problems, particularly when magnetostriction has to be measured both parallel and perpendicular to the field.

Tilted plates, with the sample placed in a hole in the lower capacitance plate, enabled the cell size to be minimised. The tilted plate construction gives good sensitivity with a sensor diameter of only  $22\text{ mm}$  and a height of  $14\text{ mm}$  (Figure 3). Careful design has provided a high dynamic measurement range, minimised the temperature dependence of the zero signal, ensured excellent thermal equilibrium, and prevented

any possibility of the capacitor plates touching. Mechanical fixing of the sample, using Cu-Be springs, enables the stress on the sample to be adjusted to a known value, and the absence of glue, for fixing or insulation, eliminates any associated capacitance drift. The use of silver for essential parts of the dilatometer halves the




Figure 1: Miniature Capacitance Dilatometer




Figure 2: Dr Martin Rotter, Prof Michael Loewenhaupt, Dr Mathias Doerr

Dr Martin Rotter and colleagues at the University of Technology, Dresden and the Technical University, Vienna, have successfully developed a very small capacitive sensor (Figure 1) for measuring thermal expansion and magnetostriction of small and irregular shaped samples[1]. The sensor can measure a range of sample lengths over a wide range of temperatures from  $0.3\text{ K}$  to  $200\text{ K}$  and in high magnetic fields up to  $15\text{ Tesla}$ . This opens up a number of possible applications, such as the investigation of crystal field effects on the magnetoelastic properties of single crystals and structural phase transitions.

#### References:

[1] M. Rotter, H. Müller, E. Gratz, M. Doerr, M. Loewenhaupt *Rev. Sci. Instr.* 692742 (1998)

#### Contact:

M. Rotter, M. Doerr, M. Loewenhaupt, Institut für Angewandte Physik und Didaktik, University of Technology, Dresden. H. Müller, E. Gratz, Institut für Experimentalphysik, Technical University, Vienna.



heat capacity compared to OFHC (oxygen free high conductivity) copper and eliminates nuclear heat capacity at low temperatures and high magnetic fields.

Although more difficult to calibrate than normal parallel plate dilatometers, Dr Rotter feels that the advantages and potential of the tilted plate sensor more than compensate for the extra effort required.




Figure 3: Arrangement of the dilatometer, temperature insert and magnet system. The full and dashed lines show the coil positions parallel and perpendicular to the magnetic field respectively.

### Oxford Instruments Superconductivity

UK  
Tubney Woods, Abingdon  
Oxfordshire OX13 5QX  
Tel: +44 (0)1865 393 311  
Fax: +44 (0)1865 393 333

Germany  
Otto-von Guericke Ring 10  
D-65205 Wiesbaden  
Tel: +49 6122 937 171  
Fax: +49 6122 937 175

Italy  
Via Leone Tolstoj 86  
20098 San Giuliano  
Milanese  
Milan  
Tel: +39 02 98 2531  
Fax: +39 02 98 24407

Japan  
Haseaman Building  
2-11-6 Tomioka, Koto-ku  
Tokyo 135-0047  
Tel: +81 03 5245 3261  
Fax: +81 03 5245 4472

Spain  
Avda Mata Piñonera, 2  
28700 San Sebastian de los Reyes  
Madrid  
Tel: +34 91 659 0740  
Fax: +34 91 654 6794

U.S.A.  
130A Baker Ave. Ext.  
Concord, MA 01742-2121  
Tel: +1 978 369 9933  
Fax: +1 978 369 6616

E-mail:  
superconductivity@oxinst.co.uk

Visit our web site at  
[www.oxford-instruments.com](http://www.oxford-instruments.com)

This publication is the copyright of Oxford Instruments Superconductivity Limited and provides outline information only which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or be regarded as a representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice, the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trade marks and registrations

© Oxford Instruments Superconductivity Limited, 2002. All rights reserved.

Reference No: OI 600/0502

UK sites of Oxford Instruments Superconductivity Limited, operate a Quality Management System approved to the requirements of BS EN ISO 9001.



Certificate No Q4118

